

Neutrophil collagenase

Catalog No: #AP85566



Package Size: #AP85566-1 50ug #AP85566-2 100ug #AP85566-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neutrophil collagenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Collagenase 2,Matrix metalloproteinase-8,MMP-8
Accession No.	Uniprot ID:O70138
Uniprot	O70138
GeneID	17394;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mmp8

Note: This product is for in vitro research use only