

## Neudesin

Catalog No: #AP85569



Package Size: #AP85569-1 50ug #AP85569-2 100ug #AP85569-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Neudesin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cell immortalization-related protein 2,Neuron-derived neurotrophic factor,Protein GIG47,Secreted protein of unknown function,SPUF protein,CIR2,SPUF
Accession No.	Uniprot ID:Q9UMX5
Uniprot	Q9UMX5
GeneID	29937;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:NENF

Note: This product is for in vitro research use only