

NADH-cytochrome b5 reductase 3

Catalog No: #AP85570



Package Size: #AP85570-1 50ug #AP85570-2 100ug #AP85570-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	NADH-cytochrome b5 reductase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	B5R,Diaphorase-1,Dia1
Accession No.	Uniprot ID:Q9DCN2
Uniprot	Q9DCN2
GeneID	109754;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cyb5r3

Note: This product is for in vitro research use only