

Nuclear receptor corepressor 2

Catalog No: #AP85571



Package Size: #AP85571-1 50ug #AP85571-2 100ug #AP85571-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor corepressor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	N-CoR2,Silencing mediator of retinoic acid and thyroid hormone receptor,SMRT,T3 receptor-associating factor,TRAC,Thyroid-, retinoic-acid-receptor-associated corepressor,Smrt
Accession No.	Uniprot ID:Q9WU42
Uniprot	Q9WU42
GeneID	20602;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ncor2

Note: This product is for in vitro research use only