

Lupus La protein

Catalog No: #AP85572



Package Size: #AP85572-1 50ug #AP85572-2 100ug #AP85572-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lupus La protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	La autoantigen,La ribonucleoprotein,Sjogren syndrome type B antigen,SS-B
Accession No.	Uniprot ID:P05455
Uniprot	P05455
GeneID	6741;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SSB

Note: This product is for in vitro research use only