

Phosphatidylcholine-sterol acyltransferase

Catalog No: #AP85574



Package Size: #AP85574-1 50ug #AP85574-2 100ug #AP85574-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphatidylcholine-sterol acyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lecithin-cholesterol acyltransferase, Phospholipid-cholesterol acyltransferase
Accession No.	Uniprot ID:P16301
Uniprot	P16301
GeneID	16816;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Lcat

Note: This product is for in vitro research use only