

Nitric oxide synthase, endothelial

Catalog No: #AP85576



Package Size: #AP85576-1 50ug #AP85576-2 100ug #AP85576-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nitric oxide synthase, endothelial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Constitutive NOS,cNOS,EC-NOS,Endothelial NOS,eNOS,NOS type III,NOSIII,Ecnos
Accession No.	Uniprot ID:P70313
Uniprot	P70313
GeneID	18127;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nos3

Note: This product is for in vitro research use only