

Myosin-4

Catalog No: #AP85577



Package Size: #AP85577-1 50ug #AP85577-2 100ug #AP85577-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myosin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Myosin heavy chain 2b,MyHC-2b,Myosin heavy chain 4,Myosin heavy chain IIb,MyHC-IIb,Myosin heavy chain, skeletal muscle, fetal
Accession No.	Uniprot ID:Q9Y623
Uniprot	Q9Y623
GeneID	4622;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MYH4

Note: This product is for in vitro research use only