

MAP kinase-activated protein kinase 3

Catalog No: #AP85578



Package Size: #AP85578-1 50ug #AP85578-2 100ug #AP85578-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAP kinase-activated protein kinase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAPK-activated protein kinase 3
Accession No.	Uniprot ID:Q3UMW7
Uniprot	Q3UMW7
GeneID	102626;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mapkapk3

Note: This product is for in vitro research use only