

## Mitogen-activated protein kinase kinase kinase MLT

Catalog No: #AP85580



Package Size: #AP85580-1 50ug #AP85580-2 100ug #AP85580-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mitogen-activated protein kinase kinase kinase MLT
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Leucine zipper- and sterile alpha motif kinase ZAK,MLK-like mitogen-activated protein triple kinase,Mixed lineage kinase-related kinase,MLK-related kinase,Sterile alpha motif- and leucine zipper-containing kinase AZK,Mltk
Accession No.	Uniprot ID:Q9ESL4
Uniprot	Q9ESL4
GeneID	65964;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Zak

Note: This product is for in vitro research use only