

Maternal embryonic leucine zipper kinase

Catalog No: #AP85582



Package Size: #AP85582-1 50ug #AP85582-2 100ug #AP85582-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Maternal embryonic leucine zipper kinase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | hMELK,Protein kinase Eg3,pEg3 kinase,Protein kinase PK38,hPK38,Tyrosine-protein kinase MELK,KIAA0175 |
| Accession No. | Uniprot ID:Q14680 |
| Uniprot | Q14680 |
| GeneID | 9833; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MELK

Note: This product is for in vitro research use only