

Methylated-DNA--protein-cysteine methyltransferase

Catalog No: #AP85584



Package Size: #AP85584-1 50ug #AP85584-2 100ug #AP85584-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Methylated-DNA--protein-cysteine methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	6-O-methylguanine-DNA methyltransferase,MGMT,O-6-methylguanine-DNA-alkyltransferase
Accession No.	Uniprot ID:P24528
Uniprot	P24528
GeneID	25332;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mgmt

Note: This product is for in vitro research use only