

Muscle, skeletal receptor tyrosine-protein kinase

Catalog No: #AP85587



Package Size: #AP85587-1 50ug #AP85587-2 100ug #AP85587-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Muscle, skeletal receptor tyrosine-protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Muscle-specific tyrosine-protein kinase receptor,MuSK,Nsk2
Accession No.	Uniprot ID:Q61006
Uniprot	Q61006
GeneID	18198;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Musk

Note: This product is for in vitro research use only