

# Protein LKAAEAR1

Catalog No: #AP85592



Package Size: #AP85592-1 50ug #AP85592-2 100ug #AP85592-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Protein LKAAEAR1                                       |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Hu   |
| Other Names        | LKAAEAR motif-containing protein 1,C20orf201           |
| Accession No.      | Uniprot ID:Q8TD35                                      |
| Uniprot            | Q8TD35   |
| GeneID             | 198437;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:LKAAEAR1

Note: This product is for in vitro research use only