

Microfibrillar-associated protein 5

Catalog No: #AP85593



Package Size: #AP85593-1 50ug #AP85593-2 100ug #AP85593-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Microfibrillar-associated protein 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MFAP-5, Microfibril-associated glycoprotein 2, MAGP-2, Magp2
Accession No.	Uniprot ID:Q9QZJ6
Uniprot	Q9QZJ6
GeneID	50530;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mfap5

Note: This product is for in vitro research use only