

## Myosin-4

Catalog No: #AP85599



Package Size: #AP85599-1 50ug #AP85599-2 100ug #AP85599-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myosin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Myosin heavy chain 2b,MyHC-2b,Myosin heavy chain 4
Accession No.	Uniprot ID:Q5SX39
Uniprot	Q5SX39
GeneID	17884;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Myh4

Note: This product is for in vitro research use only