

## Meprin A subunit alpha

Catalog No: #AP85601



Package Size: #AP85601-1 50ug #AP85601-2 100ug #AP85601-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Meprin A subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Endopeptidase-2,N-benzoyl-L-tyrosyl-P-amino-benzoic acid hydrolase subunit alpha,PABA peptide hydrolase,PPH alpha
Accession No.	Uniprot ID:Q16819
Uniprot	Q16819
GeneID	4224;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MEP1A

Note: This product is for in vitro research use only