

Mixed lineage kinase domain-like protein

Catalog No: #AP85602



Package Size: #AP85602-1 50ug #AP85602-2 100ug #AP85602-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mixed lineage kinase domain-like protein |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Accession No. | Uniprot ID:Q9D2Y4 |
| Uniprot | Q9D2Y4 |
| GeneID | 74568; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MIK1

Note: This product is for in vitro research use only