

72 kDa type IV collagenase

Catalog No: #AP85603



Package Size: #AP85603-1 50ug #AP85603-2 100ug #AP85603-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	72 kDa type IV collagenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	72 kDa gelatinase,Gelatinase A,Matrix metalloproteinase-2,MMP-2,TBE-1,CLG4A
Accession No.	Uniprot ID:P08253
Uniprot	P08253
GeneID	4313;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MMP2

Note: This product is for in vitro research use only