

Myosin-2

Catalog No: #AP85605



Package Size: #AP85605-1 50ug #AP85605-2 100ug #AP85605-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myosin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Myosin heavy chain 2,Myosin heavy chain 2a,MyHC-2a,Myosin heavy chain, skeletal muscle, adult 2
Accession No.	Uniprot ID:Q076A7
Uniprot	Q076A7
GeneID	608242;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MYH2

Note: This product is for in vitro research use only