

Canalicular multispecific organic anion transporter 2

Catalog No: #AP85606



Package Size: #AP85606-1 50ug #AP85606-2 100ug #AP85606-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Canalicular multispecific organic anion transporter 2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | ATP-binding cassette sub-family C member 3, Multidrug resistance-associated protein 3, Cmoat2, Mrp3 |
| Accession No. | Uniprot ID: B2RX12 |
| Uniprot | B2RX12 |
| GeneID | 76408; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Abcc3

Note: This product is for in vitro research use only