

DNA mismatch repair protein Msh2

Catalog No: #AP85607



Package Size: #AP85607-1 50ug #AP85607-2 100ug #AP85607-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	DNA mismatch repair protein Msh2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MutS protein homolog 2
Accession No.	Uniprot ID:P54275
Uniprot	P54275
GeneID	81709;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Msh2

Note: This product is for in vitro research use only