

Myc proto-oncogene protein

Catalog No: #AP85609



Package Size: #AP85609-1 50ug #AP85609-2 100ug #AP85609-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myc proto-oncogene protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Class E basic helix-loop-helix protein 39,bHLHe39,Proto-oncogene c-Myc,Transcription factor p64,BHLHE39
Accession No.	Uniprot ID:P01106
Uniprot	P01106
GeneID	4609;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MYC

Note: This product is for in vitro research use only