

UPF0160 protein MYG1, mitochondrial

Catalog No: #AP85610



Package Size: #AP85610-1 50ug #AP85610-2 100ug #AP85610-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	UPF0160 protein MYG1, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	C12orf10
Accession No.	Uniprot ID:Q641W2
Uniprot	Q641W2
GeneID	300258;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Myg1

Note: This product is for in vitro research use only