

## Myosin regulatory light polypeptide 9

Catalog No: #AP85611



Package Size: #AP85611-1 50ug #AP85611-2 100ug #AP85611-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myosin regulatory light polypeptide 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	20 kDa myosin light chain,LC20,Myosin regulatory light chain 2, smooth muscle isoform,Myosin regulatory light chain 9,MYRL2
Accession No.	Uniprot ID:P29269
Uniprot	P29269
GeneID	100157760;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MYL9

Note: This product is for in vitro research use only