

Norrin

Catalog No: #AP85612



Package Size: #AP85612-1 50ug #AP85612-2 100ug #AP85612-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Norrin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Norrie disease protein homolog,Ndph
Accession No.	Uniprot ID:P48744
Uniprot	P48744
GeneID	17986;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ndp

Note: This product is for in vitro research use only