

Neurofilament medium polypeptide

Catalog No: #AP85613



Package Size: #AP85613-1 50ug #AP85613-2 100ug #AP85613-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurofilament medium polypeptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	NF-M, 160 kDa neurofilament protein, Neurofilament 3, Neurofilament triplet M protein, Nef3, Nfm
Accession No.	Uniprot ID:P12839
Uniprot	P12839
GeneID	24588;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Nefm

Note: This product is for in vitro research use only