

Dual specificity mitogen-activated protein kinase kinase 2

Catalog No: #AP85614

Package Size: #AP85614-1 50ug #AP85614-2 100ug #AP85614-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dual specificity mitogen-activated protein kinase kinase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	MAP kinase kinase 2,MAPKK 2,ERK activator kinase 2,MAPK/ERK kinase 2,MEK 2,MEK2
Accession No.	Uniprot ID:Q1HG70
Uniprot	Q1HG70
GeneID	611939;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MAP2K2

Note: This product is for in vitro research use only