

## Cation-independent mannose-6-phosphate receptor

Catalog No: #AP85615



Package Size: #AP85615-1 50ug #AP85615-2 100ug #AP85615-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cation-independent mannose-6-phosphate receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CI Man-6-P receptor,300 kDa mannose 6-phosphate receptor,MPR 300,Insulin-like growth factor 2 receptor,Insulin-like growth factor II receptor,IGF-II receptor,M6P/IGF2 receptor,M6P/IGF2R,CD222,MPRI
Accession No.	Uniprot ID:P11717
Uniprot	P11717
GeneID	3482;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IGF2R

Note: This product is for in vitro research use only