

KiSS-1 receptor

Catalog No: #AP85616



Package Size: #AP85616-1 50ug #AP85616-2 100ug #AP85616-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KiSS-1 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	KiSS-1R,G-protein coupled receptor 54,G-protein coupled receptor OT7T175,hOT7T175,Hypogonadotropin-1,Kisspeptins receptor,Metastin receptor,AXOR12,GPR54
Accession No.	Uniprot ID:Q969F8
Uniprot	Q969F8
GeneID	84634;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KiSS1R

Note: This product is for in vitro research use only