

MAP kinase-activated protein kinase 3

Catalog No: #AP85619



Package Size: #AP85619-1 50ug #AP85619-2 100ug #AP85619-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAP kinase-activated protein kinase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAPK-activated protein kinase 3,Chromosome 3p kinase,3pK
Accession No.	Uniprot ID:Q16644
Uniprot	Q16644
GeneID	7867;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MAPKAPK3

Note: This product is for in vitro research use only