

Muscleblind-like protein 1

Catalog No: #AP85620



Package Size: #AP85620-1 50ug #AP85620-2 100ug #AP85620-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Muscleblind-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	RCJMB04_8n15
Accession No.	Uniprot ID:Q5ZKW9
Uniprot	Q5ZKW9
GeneID	425033;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MBNL1

Note: This product is for in vitro research use only