

Hormone-sensitive lipase

Catalog No: #AP85623



Package Size: #AP85623-1 50ug #AP85623-2 100ug #AP85623-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hormone-sensitive lipase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	HSL
Accession No.	Uniprot ID:Q68J42
Uniprot	Q68J42
GeneID	397583;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LIPE

Note: This product is for in vitro research use only