

Serine/threonine-protein kinase LMTK3

Catalog No: #AP85624



Package Size: #AP85624-1 50ug #AP85624-2 100ug #AP85624-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase LMTK3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Lemur tyrosine kinase 3,KIAA1883,TYKLM3
Accession No.	Uniprot ID:Q96Q04
Uniprot	Q96Q04
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LMTK3

Note: This product is for in vitro research use only