

MAGUK p55 subfamily member 6

Catalog No: #AP85625



Package Size: #AP85625-1 50ug #AP85625-2 100ug #AP85625-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAGUK p55 subfamily member 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Veli-associated MAGUK 1,VAM-1,VAM1
Accession No.	Uniprot ID:Q9NZW5
Uniprot	Q9NZW5
GeneID	51678;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MPP6

Note: This product is for in vitro research use only