

Membrane metallo-endoropeptidase-like 1

Catalog No: #AP85627



Package Size: #AP85627-1 50ug #AP85627-2 100ug #AP85627-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Membrane metallo-endoropeptidase-like 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Membrane metallo-endoropeptidase-like 2,NEP2(m),Nepilysin II,NEPII,Nepilysin-2,NEP2,MELL1,MMEL2,NEP2 |
| Accession No. | Uniprot ID:Q495T6 |
| Uniprot | Q495T6 |
| GeneID | 79258; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MMEL1

Note: This product is for in vitro research use only