

# 72 kDa type IV collagenase

Catalog No: #AP85629



Package Size: #AP85629-1 50ug #AP85629-2 100ug #AP85629-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	72 kDa type IV collagenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	72 kDa gelatinase,Gelatinase A,Matrix metalloproteinase-2,MMP-2
Accession No.	Uniprot ID:Q90611
Uniprot	Q90611
GeneID	386583;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MMP2

Note: This product is for in vitro research use only