

## Low-density lipoprotein receptor-related protein 2

Catalog No: #AP85632



Package Size: #AP85632-1 50ug #AP85632-2 100ug #AP85632-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Low-density lipoprotein receptor-related protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LRP-2,Glycoprotein 330,gp330,Megalin
Accession No.	Uniprot ID:A2ARV4
Uniprot	A2ARV4
GeneID	14725;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Lrp2

Note: This product is for in vitro research use only