

Myogenin

Catalog No: #AP85637



Package Size: #AP85637-1 50ug #AP85637-2 100ug #AP85637-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myogenin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Myogenic factor
Accession No.	Uniprot ID:P17920
Uniprot	P17920
GeneID	374004;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MYOG

Note: This product is for in vitro research use only