

Nuclear receptor corepressor 1

Catalog No: #AP85646



Package Size: #AP85646-1 50ug #AP85646-2 100ug #AP85646-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor corepressor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	N-CoR,KIAA1047
Accession No.	Uniprot ID:O75376
Uniprot	O75376
GeneID	9611;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NCOR1

Note: This product is for in vitro research use only