

## Neurofilament medium polypeptide

Catalog No: #AP85649



Package Size: #AP85649-1 50ug #AP85649-2 100ug #AP85649-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Neurofilament medium polypeptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	NF-M, 160 kDa neurofilament protein, Neurofilament 3, Neurofilament triplet M protein, NEF3
Accession No.	Uniprot ID:P08552
Uniprot	P08552
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: NEFM

Note: This product is for in vitro research use only