

Lysosome-associated membrane glycoprotein 2

Catalog No: #AP85651



Package Size: #AP85651-1 50ug #AP85651-2 100ug #AP85651-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysosome-associated membrane glycoprotein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LAMP-2,CD107 antigen-like family member B,LGP-96,CD107b
Accession No.	Uniprot ID:P13473
Uniprot	P13473
GeneID	3920;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LAMP2

Note: This product is for in vitro research use only