

Omega-3 fatty acid receptor 1

Catalog No: #AP85653

Package Size: #AP85653-1 50ug #AP85653-2 100ug #AP85653-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--------------------------------------------------------|
| Product Name | Omega-3 fatty acid receptor 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | G-protein coupled receptor 120,Gpr120 |
| Accession No. | Uniprot ID:Q2AC31 |
| Uniprot | Q2AC31 |
| GeneID | 294075; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:O3far1

Note: This product is for in vitro research use only