

Endonuclease 8-like 3

Catalog No: #AP85654



Package Size: #AP85654-1 50ug #AP85654-2 100ug #AP85654-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endonuclease 8-like 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DNA glycosylase FPG2,DNA glycosylase/AP lyase Neil3,Endonuclease VIII-like 3,Nei-like protein 3
Accession No.	Uniprot ID:Q8TAT5
Uniprot	Q8TAT5
GeneID	55247;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NEIL3

Note: This product is for in vitro research use only