

MAP kinase-activated protein kinase 3

Catalog No: #AP85655



Package Size: #AP85655-1 50ug #AP85655-2 100ug #AP85655-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAP kinase-activated protein kinase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MAPK-activated protein kinase 3
Accession No.	Uniprot ID:Q3SYZ2
Uniprot	Q3SYZ2
GeneID	615215;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MAPKAPK3

Note: This product is for in vitro research use only