

Nidogen-1

Catalog No: #AP85657



Package Size: #AP85657-1 50ug #AP85657-2 100ug #AP85657-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nidogen-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NID-1,Entactin,Ent
Accession No.	Uniprot ID:P10493
Uniprot	P10493
GeneID	18073;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Nid1

Note: This product is for in vitro research use only