

Glycoprotein endo-alpha-1,2-mannosidase

Catalog No: #AP85658



Package Size: #AP85658-1 50ug #AP85658-2 100ug #AP85658-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glycoprotein endo-alpha-1,2-mannosidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Endo-alpha mannosidase,Mandaselin
Accession No.	Uniprot ID:Q5SRI9
Uniprot	Q5SRI9
GeneID	79694;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MANEA

Note: This product is for in vitro research use only