

LIM domain kinase 2

Catalog No: #AP85664



Package Size: #AP85664-1 50ug #AP85664-2 100ug #AP85664-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LIM domain kinase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LIMK-2
Accession No.	Uniprot ID:O54785
Uniprot	O54785
GeneID	16886;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Limk2

Note: This product is for in vitro research use only