

Dual specificity mitogen-activated protein kinase kinase 5

Catalog No: #AP85667

Package Size: #AP85667-1 50ug #AP85667-2 100ug #AP85667-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dual specificity mitogen-activated protein kinase kinase 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAP kinase kinase 5,MAPK/ERK kinase 5,MEK 5,Mek5,Mkk5,Prkmk5
Accession No.	Uniprot ID:Q9WVS7
Uniprot	Q9WVS7
GeneID	23938;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Map2k5

Note: This product is for in vitro research use only