

# Macrophage mannose receptor 1

Catalog No: #AP85668

Package Size: #AP85668-1 50ug #AP85668-2 100ug #AP85668-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Macrophage mannose receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MMR,C-type lectin domain family 13 member D,C-type lectin domain family 13 member D-like,Human mannose receptor,hMR,Macrophage mannose receptor 1-like protein 1,CD206,CLEC13D,CLEC13DL,MRC1L1
Accession No.	Uniprot ID:P22897
Uniprot	P22897
GeneID	4360;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MRC1

Note: This product is for in vitro research use only